

FHIRDBPKG FoundationProduction Configuration

- [Configuration Summary](#)
- [Business Hosts](#)

Configuration Summary

This report lists the configured items in the production and the values of their most important configuration settings.

Platform

The following table describes the instance that generated this production report.

Item	Value
Report Time and Date	2022-09-30 12:47:32
Version	IRIS for UNIX (Red Hat Enterprise Linux 7 for x86-64) 2022.1 (Build 209U)

Production

The following settings apply to the production as a whole.

Setting	Value
ActorPoolSize	1
Description	
LogGeneralTraceEvents	false
Name	FHIRDBPKG.FoundationProduction

Interfaces

Business hosts in the production specify their interactions with one another as follows. This table provides a summary only. For additional details about specific items in this production, refer to Business Hosts and other sections in this report.

Business Host	Setting	Value
FromDropbox	TargetConfigNames	MFT.FHIRProcess, ToBox
FromBox	TargetConfigNames	MFT.FHIRProcess, ToDropbox

External Systems

Business hosts in the production connect with external systems as follows.

Business Host	Setting	Value
FromDropbox	DeleteFromServer	1
FromBox	DeleteFromServer	2

Business Hosts

This section describes the configuration settings for business hosts in the production. Settings that can be true or false also accept the values 1 or 0, respectively.

FromBox

The following tables list configuration settings that apply to the business host FromBox.

Setting	Value
Category	Box, MFT
ClassName	EnsLib.MFT.Service.Passthrough
Comment	Ingest FHIR from Dropbox
Enabled	true

Foreground	false
LogTraceEvents	false
Name	FromBox
PoolSize	1
Schedule	
DeleteFromServer	2
FilesToRetrieve	
IncludeSubFolders	1
MFTConnectionName	FHIRBox
MFTSourceFolders	/FHIRBox
TargetConfigNames	MFT.FHIRProcess, ToDropbox

FromDropbox

The following tables list configuration settings that apply to the business host FromDropbox.

Setting	Value
Category	Dropbox, MFT
ClassName	EnsLib.MFT.Service.Passthrough
Comment	Ingest FHIR from Dropbox
Enabled	true
Foreground	false
LogTraceEvents	false
Name	FromDropbox
PoolSize	1
Schedule	
DeleteFromServer	1
MFTConnectionName	FHIRDrop
MFTSourceFolders	/InFHIR
TargetConfigNames	MFT.FHIRProcess, ToBox

MFT.FHIRProcess

The following tables list configuration settings that apply to the business host MFT.FHIRProcess.

Setting	Value
Category	MFT, Box, Dropbox
ClassName	ZMSP.MFT.FHIRProcess
Comment	Process FHIR
Enabled	true
Foreground	false
LogTraceEvents	false
Name	MFT.FHIRProcess
PoolSize	1
Schedule	

ToBox

The following tables list configuration settings that apply to the business host ToBox.

Setting	Value
Category	Dropbox, Box, MFT
ClassName	EnsLib.MFT.Operation.Passthrough
Comment	
Enabled	true
Foreground	false

LogTraceEvents	false
Name	ToBox
PoolSize	1
Schedule	
DefaultMFTFolder	/FromDropbox
MFTConnectionName	FHIRBox

ToDropbox

The following tables list configuration settings that apply to the business host ToDropbox.

Setting	Value
Category	Box, Dropbox, MFT
ClassName	EnsLib.MFT.Operation.Passthrough
Comment	
Enabled	true
Foreground	false
LogTraceEvents	false
Name	ToDropbox
PoolSize	1
Schedule	
DefaultMFTFolder	/FromBox
MFTConnectionName	FHIRDrop