

Is there an API for determining the size of a routine?

Published on InterSystems Developer Community (<https://community.intersystems.com>)

Question

[Renee Cannon](#) · Mar 28, 2016

Is there an API for determining the size of a routine?

Is there an API I can call to determine the size of a routine in a CACHE.DAT file so that I don't have to directly reference the ^ROUTINE global?

[#Databases](#) [#Object Data Model](#) [#ObjectScript](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/there-api-determining-size-routine>