
Question

[Ricardo Paiva](#) · Jul 13, 2016

Overriding soap body

We can override soap body by using WriteSOAPBodyMethod. Something like:

```
Method Online(ByRef Obj1 As myClass1, Obj2 As myClass2, Output ObjOutput1 As myClass1Output(XMLNAME="Obj")) [ ..., WebMethod ]
```

```
{  
  
    ...  
  
    s ..WriteSOAPBodyMethod="override"  
  
    d (..WebMethod("Online","ProcessRequest")).Invoke($this,"Online",.Obj1,.Obj2,. ObjOutput1)  
  
    s ..WriteSOAPBodyMethod=""  
  
    ...  
  
}
```

```
Method override(proxy As %SOAP.ProxyDescriptor, tag As %String)
```

```
{  
  
    ...  
  
}
```

In order to generate soap body, how can one access Obj1 and Obj2 from within the override method?

[#SOAP #Caché](#)

Source URL: <https://community.intersystems.com/post/overriding-soap-body>