

---

## Article

[Mike Kadow](#) · Sep 2, 2016 1m read

## Quickest way to backup routines

I have found this one line command from the Terminal is the quickest way to backup my routines.

In the Namespace the routines reside:

```
w $SYSTEM.OBJ.Export("MYRTNS*.MAC","C:/Users/Username/Desktop/MYRTNS.BCK-9-2-2016-8AM")
```

[#Caché](#) [#Terminal](#) [#Backup](#)

---

Source URL: <https://community.intersystems.com/post/quickest-way-backup-routines>