

---

Question

[Mikhail Khomenko](#) · Oct 30, 2016

## ISCAgent debug

Good day!

Is there an opportunity to debug the ISCAgent behaviour (in Linux)? ISCAgent uses `/etc/iscagent/iscagent.conf` as configuration file, but in docs I've found a description only for two parameters (port number and interface - <http://docs.intersystems.com/latest/csp/docbook/DocBook.UI.Page.cls?KEY=...>). Other parameters I can see by running `/usr/local/etc/cachesys/ISCAgentCtrl status`:

```
application_server.interface_address=*
application_server.port=2188
daemonized=yes
mirroring=on
pid=41850
running=yes
version=2015.2.0.664.3
```

And in `/var/log` there is a file called `iscagentconsole.log`, which is always empty.

So, the questions are the next. How can I debug ISCAgent behaviour excluding `strace` usage? Can I add debug level parameter to `iscagent.conf` and see debug lines in `iscagentconsole.log`?

Thanks a lot!

[#System Administration](#) [#Mirroring](#) [#Caché](#)

---

Source URL: <https://community.intersystems.com/post/iscagent-debug>