

Question

[Bapu Hirave](#) · Nov 14, 2016

Soap Response performance

I have a soap webservice extended from `%SOAP.WebService` class

Lets say i have web method and have return type of following classes (either one of them)

Class `ReturnTypeA`

```
{  
    Property SelectedDoses As %XML.DataSet;  
}
```

Class `ReturnTypeB`

```
{  
    Property SelectedDoses As list Of WebServices.DataTypes.Dose;  
}
```

Which return type class has better performance while in terms of xml serialization ? A or B

my client application is in C#

[#Object Data Model](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/soap-response-performance>