

Question

[Kurro Lopez](#) · Jan 26, 2017

SOAP.OutboundAdapter HttpResponseMessage doesn't handle large response [SOLVED]

Hi all,

Sorry to bother 

I have a Business Object with SOAP.OutboundAdapter and it gets the response from a WebApi.

Following the sample [Creating REST Services and Clients with Ensemble - Developing a REST Operation](#), I've done the call and it retrieves a JSon with the information. The result is a JSon with a Base64 content (at least 15000 bytes)

I've put a trace to check the content with \$\$\$TRACE(tHttpResponse.Data.Read() I have a cut result (to 1200 character), so the JSonStreamToObject is not able to allocate the full result into a object.

is there any way to handle more character in HttpResponseMessage or is about SOAP.OutboundAdapter configuration?

Best regards

[#Caché](#) [#Ensemble](#) [#Business Operation](#)

Source

URL: <https://community.intersystems.com/post/soapoutboundadapter-httpresponse-doesnt-handle-large-response-solved>