
Question

[Steve Pisani](#) · Mar 31, 2017

Programatic User Authentication against windows AD

Hi,

Does calling the BIND method of %SYS.LDAP, with the username, domain and password of the user that needs to be authenticated- the right way to authenticate him/her ?

Also - am I correct in assuming that something like this is independant to (and I don't need to specify setting for), System Security -> LDAP Options

Thanks

Steve

[#Microsoft Windows](#) [#Authentication](#) [#LDAP](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/programatic-user-authentication-against-windows-ad>