
Question

[Yaniv Ben Malka](#) · Oct 10, 2017

Index Recommended Approaches

Hi,

I have a class with around 400k lines and 60 columns. Class storage is Cache SQL storage (Mapped from a global).

I want to create multiple indices on certain fields.

I am familiar with two approaches:

1. Create a new map (Index type) on a pointer global.
2. Create a bitmap index

Which approach is more recommended to be used in the case I described? If there are any other approaches, I will be happy to hear.

Thanks :)

[#Globals](#) [#Indexing](#) [#SQL](#) [#Cache](#)

Source URL: <https://community.intersystems.com/post/index-recommended-approaches>