
Question

[Maks Atygaev](#) · Nov 27, 2017

\$lb vs \$c for byte array

Dear community!

I'm currently working on MessagePack implementation by using Cache ObjectScript.

I try to choose suitable representation for binary array.

Should I use \$lb or \$c or smth else for it?

<https://msgpack.org/index.html>

<https://github.com/msgpack/msgpack/blob/master/spec.md>

[#ObjectScript](#) [#Cache](#)

Source URL: <https://community.intersystems.com/post/lb-vs-c-byte-array>