

Question

[Sébastien Demoustiez](#) · Dec 14, 2017

Get a custom response server header with Cache

Hello,


I'm using Cache 2016.

I try to get a property set in my server response header.

I tried to get it in %request and %response but without success ...

Some one have an idea ?

Here an example of what i try to get :



The screenshot shows the 'Headers' tab of a web browser's developer tools. The 'Response Headers' section is highlighted with a red box. The headers listed are:

- Cache-Control: no-cache
- Connection: Keep-Alive
- Content-Encoding: gzip
- Content-Length: 852
- Content-Security-Policy: default-src 'self' https://ldp.iamfas.belgium.be https://www.youtube.com http://intranetappadm; script-src 'self' 'unsafe-inline' 'unsafe-inline'
- Content-Type: text/html; charset=utf-8
- Date: Thu, 14 Dec 2017 13:23:53 GMT
- Expires: Thu, 29 Oct 1998 17:04:19 GMT
- Keep-Alive: timeout=5, max=298
- Pragma: no-cache
- Referrer-Policy: origin
- Set-Cookie: CSPSESSIONID=SP-80-UP=-0000010100003ryK2NHA4HT8u9RTck_4X1jFhCvpJQk20nJwm; path=/; httpOnly
- Set-Cookie: CSPSERVERID=A125Xyi3; path=/; httpOnly
- X-Content-Type-Options: nosniff
- X-Frame-Options: SAMEORIGIN
- X-Permitted-Cross-Domain-Policies: master-only
- X-XSS-Protection: 1; mode=block

The 'Request Headers' section is also visible below, showing:

- Accept: text/html,application/xhtml+xml,application/xml;q=0.9,image/webp,image/apng,*/*;q=0.8
- Accept-Encoding: gzip, deflate, br
- Accept-Language: fr-FR, fr;q=0.9,en-US;q=0.8,en;q=0.7,nl;q=0.6
- Cache-Control: max-age=0
- Connection: keep-alive
- Cookie: CSPSESSIONID=SP-80-UP=-0000010100003ryK2NHA4HH5bNmUw5JkdTRpC_97gavjrpZ2zR; CSPSERVERID=A125Xyi3
- Host: api.hap.be
- Upgrade-Insecure-Requests: 1
- User-Agent: Mozilla/5.0 (Windows NT 6.1; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/62.0.3202.94 Safari/537.36

Thank you

Sébastien

[#Object Data Model](#) [#Frontend](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/get-custom-response-server-header-cache>