
Question

[Arun Kumar](#) · Feb 1, 2018

Performance in Cache ObjectScript

Hi Guys,

Can you please advise on the below queries.

Query 1:

Example 1:

```
S a="345",b="arun",c="kumar",d="hi",e="yello",f="orange"
```

Example 2:

```
S a="345"
```

```
S b="arun"
```

```
S c="kumar"
```

```
S d="hi"
```

```
S e="yello"
```

```
S f="orange"
```

Can you please advise me, which one is performance wise is better.

Query 2:

Example 1:

```
S:a=1 R="Arun"
```

Example 2:

```
I a=2 R="Arun"
```

Please advise me, which one is giving better performance in this.

Any lead would be appreciated.

Thanks,

Arun Kumar Durairaj.

[#Cache](#) [#ObjectScript](#) [#Performance](#)

Source URL: <https://community.intersystems.com/post/performance-cache-objectscript>