
Question

[George Hodder](#) · Feb 2, 2018

OAuth2.0 with ZAUTHENTICATE routine and Rest service

Hi.. I have an issue where we are using OAuth2.0 with the ZAUTHENTICATE routine. Once our token is validated we are using a users lan id (passed on the ID token) to find a software defined username in a Cache Global.

That is all working fine in ZAUTHENTICATE.. I am setting the software defined username in the Properties("Comment") array and wanting to reference it in the Rest Service Dispatch class.

I am finding that Properties array is not resident in the dispatch class. It is my understanding that Properties is passed to ZAUTHENTICATE by reference, but how can I invoke the dispatch class in a way that includes the array? Is there setup that needs added in the portal for web application security?

Additionally I see that ZAUTHENTICATE quits with 1/0 for success/fail.. is there a way to return more specific messaging regarding the failure to the calling web application?

Thanks,

George

[#Beginner](#) [#OAuth2](#) [#REST API](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/oauth20-zauthenticate-routine-and-rest-service>