

Question

[Guillaume Lepretre](#) · May 14, 2018

Get request header (REST API)

Hello,

In all web services, i need to my get login and token. So with Postman, i tried to call a HTTP request where I put the login/token in the header :

Authorization	Headers (3)	Body	Pre-request Script	Tests
	Key			Value
<input checked="" type="checkbox"/>	Content-Type			application/json
<input checked="" type="checkbox"/>	Login			testUtilisateur
<input checked="" type="checkbox"/>	Token			token
	New key			Value

I tried to get data from Http request header. The REST API use %CSP.REST. I tried something like that :

```
+ y %  
set login = %request.GetHeader("LOGIN")  
set login = %request.Login  
|
```

But it didn't work..

Someone can give me some example or other method ?

Regards,

[#API](#) [#Ensemble](#) [#REST API](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/get-request-header-rest-api>