

Question

[Chandra Bandi](#) · Sep 24, 2018

How to create a restful api in Intersystems cache with %CSP.REST

Hello Guys,

Can you please guide me to create a RESTful API service in our cache (using intersystems cache kit) with CSP (Cache Server Pages) object script.

```
Class REST.NativeDispatcher Extends %CSP.REST
{
XData UrlMap [XMLNamespace = "http://www.native.rest.com/urlmap" ]
{
<Routes>
  <!-- <Route Url="/:name" Method="GET" Call="displaySystem" Cors="false" /> -->
  <Route Url="/:ostype" Method="GET" Call="externalFreeze" Cors="false" />
</Routes>
}

ClassMethod externalFreeze(ostype as %String) as %Status{
    set status = ##Class(Backup.General).ExternalFreeze()
    WRITE "EXTERNAL FREEZE RESPONSE FOUND " _status
    QUIT $$$OK
}
}
```

I had written above code snippet but, I need to explore more using %CSP.REST object script language.

Can anyone suggest/ share sample code snippet to create a REST service web application using %CSP.REST object script language.

Thanks,

Chandrasekhar Reddy,

Email: cbandi@purestorage.com | +91 9880318877

[#Object Data Model](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/how-create-restful-api-intersystems-cache-csprest>