

Terminal can not be connected.

Published on InterSystems Developer Community (<https://community.intersystems.com>)

Question

[akio tanaka](#) · Sep 25, 2018

Terminal can not be connected.

Setting up CACHE preferred connection server.

server name:localhost

ip adress:127.0.0.1

port:1972

Telnetport:23

web server port:57772

web server ip adress:

CSP server instance:

comment:

Authentication method=>password

hosts

localhost 127.0.0.1

os

windows10 (virtual box)

>ifconfig

windows IP configuration

Ethernet adapter Ethernet:

Connection-specific DNS suffix:

Link local IPv6 address: fe80 :: b943: 6646: 21a8: 1a5b% 2

IPv4 address: 10.0.2.15

Subnet mask: 255.255.255.0

Default gateway: 10.0.2.2

>ping 127.0.0.1

Number of packets: transmission = 4, reception = 4, loss = 0 (0% loss)

>netstat -an

127.0.0.1:1900 :

127.0.0.1:55152 :

127.0.0.1:57967 :

nets winhttp>show proxy

Current WinHTTP proxy settings:

Direct access(no proxy server)

If you enter the preferred connection server as above and click on the cube and make it to remote system access terminal localhost,

TCP connect () failed - exception satisfied select (), Reason: (10061, 0x274 d) The connection could not be made because it was rejected by the target computer.

It will be displayed.

Is the way of setting wrong?

It would be greatly appreciated if you could tell us what you can fix and you will be able to use the terminal.

Terminal can not be connected.

Published on InterSystems Developer Community (<https://community.intersystems.com>)

Thanking you in advance.

[#Caché](#)

Source URL: <https://community.intersystems.com/post/terminal-can-not-be-connected>