

How to mimic selenium in cache object Script

Question

[Thembelani Mlalazi](#) · Dec 11, 2018

How to mimic selenium in cache object Script

Is it possible to mimic what selenium does like navigating to a site and logging in and filling out a form then logout in COS. I am trying to do that in COS using %Net.HttpRequest class or should I be using a different class the idea is to be able to call a web app login into it and fill out form and log out

[#API](#) [#Beginner](#) [#Object Data Model](#) [#ObjectScript](#) [#REST API](#) [#Testing](#) [#Cache](#)

00 2 0 7 347

[Log in or sign up to continue](#)

[Add reply](#)

Source URL: <https://community.intersystems.com/post/how-mimic-selenium-cache-object-script>