
Question

[Krishnamuthu Ve...](#) · Apr 9, 2019

What is the preferred / best way of logging the Request and Response of REST calls in Intersystems Cache ?

I want to log the input request and response to the REST call.

There is an option for SOAP web service `^ISCSOAP("Log")`.

Is there anything similar for REST services in Cache ? Or what is the preferred way of logging the REST calls.

[#Caché](#)

Source

URL:<https://community.intersystems.com/post/what-preferred-best-way-logging-request-and-response-rest-calls-intersystems-cache>