

---

Question

[William Proctor](#) · May 2, 2019

## replacement for Putty connecting to cache

Thanks for all replies in advance. We have a security vulnerability that we have to get rid of. We use Putty software to connect to cache as a terminal allowing several users to do maintenance work in cache. this uses telnet Plain text. I know that we can configure telnet to be encrypted using the super server service and I'm looking for software that can work like Putty as a terminal using encryption compatible with cache telnet encryption. If I have cache installed on my PC and setup a connection to the server using Kerberos with encryption and use the terminal option to connect to the server will it be encrypted or will the terminal telnet session still be plain text.

Authentication Method

- ☐ Password
- ☒ Kerberos

Connection Security Level

- ☐ Kerberos
- ☐ Kerberos with Packet Integrity
- ☒ Kerberos with Encryption

☒ Windows Caché Telnet Server

[#Databases](#) [#Encryption](#) [#FTP](#) [#Caché](#)

---

Source URL: <https://community.intersystems.com/post/replacement-putty-connecting-cache>