
Question

[Stas Rabkin](#) · May 19, 2019

Using VSTS As Version Control for Caché code

Hi,

recently we started to work with Atelier & VSTS in our organization .

we are looking for a way to use the VSTS as Our Version/Source control system .

I have read the article about "[Continuous Delivery of your InterSystems solution using GitLab](#)" & now trying to find a way to use a similar methodology with the VSTS that we have in our organization.

It Appears that it is possible to Create a Git Repository in the VSTS but i'm not sure how are we to use it from the atelier.

The main questions are:

1. What Is the best practice in saving cache code in VSTS GIT?(cls files?, legacy XML?...)
2. What is the right flow of importing/exporting code in/from the development repository during the development process?
3. What is the right flow of building/releasing cache code to higher environment?(dev-->Test, Test-->Prod)

we are using :

eclipse Photon Release (4.8.0)

Cache for Windows (x86-64) 2017.1.1 (Build 111U) Wed May 17 2017 15:34:56 EDT

Thanks ,

Stas Rabkin.

[#Git](#) [#Caché](#) [#Other](#)

Source URL:<https://community.intersystems.com/post/using-vsts-version-control-cach%C3%A9-code>