

---

### Question

[Robert Cemper](#) · Sep 12, 2019

## Features of \*ENSTEMP

ENSTEMP was addressed already some time ago;

[Article: Preventing Globals From Getting Journalled \(Continued from How do I Minimize My Journals\)](#)

and is also in public documentation

[Where InterSystems IRIS Stores Temporary Production Data](#)

Though I'd like to understand if this is just naming a non-journalled DB  
or are there also typical features of CACHETEMP / IRISTEMP related to it:

- automatic clean-up at system startup
- keeping Global Buffers in memory as long as possible
- late writing to the storage file

In other words is it a 2nd, 3rd, .. fully feature blown IRISTEMP ?

[#Other](#)

---

Source URL: <https://community.intersystems.com/post/features-enstemp>