
Question

[Yaron Munz](#) · Oct 16, 2019

Ensemble dynamic archive path

Hello,

Is there an elegant way to define a dynamic "Archive path" for `EnsLib.File.InboundAdapter` ?

What i mean is to have a sub-folders structure (under the root archive path) with a name similar to: `yyyymmdd`
This way the adapter "could" store the files on those sub-folders (and create a new one for the 1st file for a specific date).

When a production is processing hundreds of thousands of files every day, we might end up with an archive with huge number of files (many millions) which is then hard to handle (read, search) especially on Windows O/S.

Of course that there is always a "programmatic" way to move/distribute the files from archive root (as a background task), but before going to that direction, I would like to know, if there is an "out of the box" solution.

Regards.

Yaron

[#Ensemble](#)

Source URL: <https://community.intersystems.com/post/ensemble-dynamic-archive-path>