
Question

[Kevin Kindschuh](#) · Jan 15, 2020

How can I determine the IP address from which data is come from a TCP Inbound Adapter?

We want to be able to dispatch based on the IP address who connected to a TCP Inbound adapter.

The address is in `..Adapter.%outIPAddr` when the connection is made and reported with `$$$LOGINFIO()` but by the time `OnProcessInput()` is called the value is blank.

[#Ensemble](#)

Source

URL: <https://community.intersystems.com/post/how-can-i-determine-ip-address-which-data-come-tcp-inbound-adapter>