

Announcement

[Anastasia Dyubaylo](#) · May 13, 2020

[Open Exchange](#)

New Coding Talk: Using the IRIS Native API GitHub Template

Hi Developers,

Please welcome the new video specially recorded by InterSystems Product Manager [@Bob Kuszewski](#) for the [3rd InterSystems Online Programming Contest](#) on InterSystems IRIS Native API:

[Using the IRIS Native API GitHub Template](#)

What's the contest about?

In the last year, InterSystems added and expanded [Native APIs](#) for Java, .NET, Python, and Node.js. The Native API provides high speed access to globals, the ability to instantiate and use ObjectScript objects, and fine-grained control over transactions.

In this contest, we would like to see what kinds of interesting applications can be built using the Native API.

[Native API Contest Template](#)

To hope you get started with the contest, we've provided this template in GitHub to get you started. Just follow the instructions included there to clone the repository, build a docker container, and run your application.

We hope that you'll enjoy the new APIs and have fun hacking!

Stay tuned!

[#NET #API #C++ #Contest #Docker #GitHub #Java #Node.js #ObjectScript #Python #Video #InterSystems IRIS #InterSystems IRIS for Health #Open Exchange](#)
[Check the related application on InterSystems Open Exchange](#)

Source URL: <https://community.intersystems.com/post/new-coding-talk-using-iris-native-api-github-template>