
Announcement

[Anastasia Dyubaylo](#) · Jun 5, 2020

[Open Exchange](#)

New Coding Talk: OData Server for InterSystems IRIS

Hey Developers,

New "Coding Talk" video was specially recorded by [@Yuri Marx](#) for the [InterSystems IRIS Native API Programming Contest](#):

[OData Server for InterSystems IRIS](#)

In this video, you'll learn about the OData Server for InterSystems IRIS SQL Repositories. This project allows you to virtualize Persistent classes from InterSystems IRIS to REST OData Services.

[OData Server for IRIS](#) on Open Exchange

For any questions, please write to [@Yuri Marx](#) at yurimarx@gmail.com or in the Direct Messages on DC.

Enjoy watching this promo video!

And don't forget to vote for this project in the [IRIS Programming Contest](#)!

[#Video #Open Exchange](#)

[Check the related application on InterSystems Open Exchange](#)

Source URL: <https://community.intersystems.com/post/new-coding-talk-odata-server-intersystems-iris>