


Article

[Yuri Marx](#) · Jun 22, 2020 1m read

## Create cognitive digital service with InterSystems IRIS



Powered by  
**bizagi**  
Modeler

The InterSystems IRIS is a great platform to develop, run and consume Data Science services. IRIS can ingest data from any type, format, protocol and time using adapters. These datasets can be prepared with BPL, DTL and Object Script and stored as SQL or NoSQL data. Finally, it can be consumed by open ML algorithms inside IRIS and visualized in the IRIS dashboards. See more in:

<https://docs.intersystems.com/irislatest/csp/docbook/Doc.View.cls?KEY=PAGEdatascience>.

[#InterSystems Business Solutions and Architectures](#) [#Machine Learning \(ML\)](#) [#InterSystems IRIS](#)

Source URL: <https://community.intersystems.com/post/create-cognitive-digital-service-intersystems-iris>