

Announcement

[Anastasia Dyubaylo](#) · Jul 13, 2020

Vote for the best app in the InterSystems IRIS AI Programming Contest!

Hey Developers!

This week is a voting week for the [InterSystems IRIS AI Programming Contest!](#)

So, it's time to give your vote to the best AI- and ML-enabled solution on InterSystems IRIS!

You decide. [VOTING IS HERE](#)



How to vote? This is easy: you will have one vote, and your vote goes either in Experts Nomination or in Community Nomination.

Experts Nomination

If you are InterSystems Product Manager, or DC moderator, or Global Master from Specialist level and above cast your vote in the Expert nomination.

Community Nomination

If you ever contributed to DC (posts or replies) and this was not considered as spam you are able to vote for the applications in the Community nomination.

Voting

Registered applications

Experts

Community

identity_cache

ASP.NET Identity Cache Provider

Author: [Evgeny Shvarov](#)

Experts 0

Community 0

Vote



DeepSee Audit

Example of Deepsee usage for Caché© Audit

Author: [Semen Makarov](#)

Experts 0

Community 0

Vote



Voting takes place on the [Open Exchange Contest Page](#) and you need to sign in to [Open Exchange](#) — you can do it with your DC account credentials.

If you changed your mind, cancel the choice and give your vote to another application - you have 7 days to choose!

Contest participants are allowed to fix the bugs and make improvements to their applications during the voting week, so don't miss and subscribe to application releases!

Also, please check out the [Judgment and Voting Rules for the Contest here](#).

[#Artificial Intelligence \(AI\)](#) [#Contest](#) [#IntegratedML](#) [#Machine Learning \(ML\)](#) [#InterSystems IRIS](#) [#IRIS contest](#) [#Open Exchange](#)

Source URL: <https://community.intersystems.com/post/vote-best-app-intersystems-iris-ai-programming-contest>