
Announcement

[Anastasia Dyubaylo](#) · Jul 18, 2020

[Open Exchange](#)

New Coding Talk: Monitoring your ML model with IntegratedML

Hi Developers,

Please welcome the new video, specially recorded by [@José Pereira](#) for the [InterSystems IRIS AI Programming Contest](#):

[Monitoring your ML model with IntegratedML](#)

This video shows you a technology example of a machine learning model being monitored by a custom system monitor. E-mail alerts will be triggered when model accuracy is less than 90%.

The application is served as a docker container, so you can easily run it on your system. Code and instructions are here:

[iris-integratedml-monitor-example](#) on Open Exchange

Why is COVID-19 also dangerous for Machine Learning?

- [Part 1](#)
- [Part 2](#)

For any questions, please write to [@José Pereira](#) at joseroberto.pereira@shift.com.br or in the Direct Messages on DC.

Support this project with your vote in the [IRIS AI Programming Contest](#)!

Enjoy!

[#AI](#) [#Contest](#) [#IntegratedML](#) [#Machine Learning](#) [#Python](#) [#Video](#) [#InterSystems IRIS](#) [#Open Exchange](#)
[Check the related application on InterSystems Open Exchange](#)

Source URL: <https://community.intersystems.com/post/new-coding-talk-monitoring-your-ml-model-integratedml>