

Unable to set DG1 segment after cloning

Question



[ED Coder](#) · Sep 17, 2020

Unable to set DG1 segment after cloning

My segment is not getting set with the new value though It goes into the if statement. Am I doing something wrong?

```
ClassMethod Clone(pRequest As EnsLib.HL7.Message) As EnsLib.HL7.Message [ Final ]
```

```
s cREQ = pRequest.%ConstructClone()
```

```
s DG1 = cREQ.FindSegment("DG1",.ind)
```

```
While (ind2 != "")
```

```
{  
s diag = DG1.GetValueAt(4.1)
```

```
if diag = "0"
```

```
{  
s tmp = "No code"  
s rc=DG1.SetValueAt(tmp,4.1)  
}
```

```
}
```

```
s DG1 = cREQ.FindSegment("DG1",.ind2)
```

```
s request = cREQ  
Quit request
```

```
}
```

00 0 6 77 [3](#)

Source URL: <https://community.intersystems.com/post/unable-set-dg1-segment-after-cloning>