

Configure File Service to retrieve a file every 10 minutes

Question

[Mikael Toivonen](#) · Oct 12, 2020

Configure File Service to retrieve a file every 10 minutes

I have a file passthrough service that should check for new new files every 10 minutes, process them and then wait again for 10 minutes, even if there are new files waiting in the source directory. The "CallInterval" value feels like a wrong tool for this need since it would keep the service "locked" for a maximum of 10 minutes before releasing resources. Is there a better way than CallInterval to schedule the service to poll for files every 10 minutes?

[#Business Service](#) [#InterSystems IRIS](#)

00 2 0 2 92

[Log in or sign up to continue](#)
[Add reply](#)

Source URL: <https://community.intersystems.com/post/configure-file-service-retrieve-file-every-10-minutes>