
Question

[Scott Fadden](#) · Oct 27, 2020

Configure IRIS using Python

I am trying to configure IRIS using python. The first task I wanted to do it change the journal directory. I cannot figure out the right way to work with the Config.Journal class with python. I can get the text to query the current setting.

```
returnValue = dbnative.classMethodValue("Config.Journal","GetList")  
print(returnValue )
```

But I don't see a straightforward way to parse the \$List returned by GetList() I was hoping there was an easier way.

How do I

1. Query configuration parameters?
2. Set configuration parameters?

Using python?

[#Python](#) [#InterSystems IRIS for Health](#)

Source URL: <https://community.intersystems.com/post/configure-iris-using-python>