

---

Question

[Phung Thoa](#) · Nov 2, 2020

## Where to save the local cache code after connecting the IRIS server using Visual Studio Code

Hi all,

By using VS code, I want to copy multiple files to local directory, then sync them to server, but I can't find where local directory is stored?

In Atelier IDE, I could open a view of server and can copy the code to the project, then code is saved in the local directory. When I copy multiple files to the local directory -> choose sync in Atelier -> the code will be synchronized to the server

How to do this from VS Code?

Thanks!

[#Caché](#) [#InterSystems IRIS](#) [#VSCode](#)

---

Source

URL:<https://community.intersystems.com/post/where-save-local-cache-code-after-connecting-iris-server-using-visual-studio-code>