

---

Question

[Cindy Fawcett](#) · Apr 12, 2021

## API Manager - transform header from basic auth to bearer token

I am working with InterSystems API Manager product and I'm trying to solve an authentication problem. In our implementation, API Manager accepts basic authentication. One of the endpoints accepts OAuth2. The request coming into API Manager is FHIR and needs to have basic auth in the header. The OAuth token is passed to API Manager as a parameter. I want to send FHIR to the endpoint with a Bearer Token in the header instead of the basic auth that came into API Manager. Has anyone used a transformation plug-in to solve a similar issue? I'm trying to solve this problem in API Manager rather than the API endpoint.

[#InterSystems API Manager \(IAM\)](#) [#InterSystems IRIS for Health](#)

---

Source URL: <https://community.intersystems.com/post/api-manager-transform-header-basic-auth-bearer-token>