

---

Question

[Thomas Schroyen](#) · Jul 14, 2021

## creation of a bidirectionnial connector

Hello everyone,

I have a question and if you have some code, it would be lovely.

I had to create a connector which is biredictionnal. What does I mean by this?

I have to send ADT to an certain ip +port by TCP but I have to receive query A19 by this same ip+port.

Is it possible to create such a connector that can do both?

I read about pipe but i think it's too low-leveled .

Thanks,

Thomas

[#HealthShare](#)

Product version: HealthShare 2018.1

---

Source URL: <https://community.intersystems.com/post/creation-bidirectionnial-connector>