
Question

[Evgeny Shvarov](#) · Nov 10, 2021

How to Implement Bearer Authentication with InterSystems IRIS REST API?

Hi folks!

Could you please share your experience on how do you create REST API with InterSystems IRIS that uses bearer authentication?

How do you generate tokens? How do you maintain it (how much time tokens exist?).

Thanks in advance!

[#Authentication](#) [#REST API](#) [#InterSystems IRIS](#)

Product version: IRIS 2021.1

Source

URL: <https://community.intersystems.com/post/how-implement-bearer-authentication-intersystems-iris-rest-api>