
Question

[phillip jack](#) · Feb 2, 2022

Saving a serial object using Result set

Hi Folks,

I created a persistent class as below

```
Class myclass.DataBase Extends %Persistent
{
Property ID As %String;

Property Skill As list Of myclass.SerialTablelist;
}
```

and Created another Serial class as

```
Class myclass.SerialTablelist Extends %SerialObject
{
Property PSkill As %String;

Property OSkill As %String;
}
```

Now I will save the id as below

```
do rs.Prepare("Insert into myclass.DataBase(ID)VALUES(?)")
do rs.Execute(ID)
```

I am facing issue with the AS LIST OF

Could any one guide me how to save the Skill (PSkill ,OSkill) using resultset (I am able to)

[#Databases](#) [#InterSystems Business Solutions and Architectures](#) [#ObjectScript](#) [#Caché](#)

Source URL: <https://community.intersystems.com/post/saving-serial-object-using-result-set>