
Discussion

[Vikram Rajput](#) · Feb 15, 2022

IRIS Native API for .NET performance issue

I am surprised with the performance (poor) of Native API for .NET

Following code (where ^testGlobal,node1 has 50 records) takes 5 seconds to run. Anyone else think there is room for the improvement?

```
string global = "^testGlobal";
object[] Subs = new object[1];
Subs[0] = "node1";

try
{
    IRISIterator iter = iris.GetIRISIterator(global, Subs);
    foreach (var item in iter)
    {
        Console.WriteLine((string)(iter.CurrentSubscript));
    }
}
catch (Exception ex)
{
}
}
```

[#.NET](#) [#Performance](#) [#InterSystems](#) [IRIS](#)

Source URL: <https://community.intersystems.com/post/iris-native-api-net-performance-issue>