
Question

[Carl Emberger](#) · Feb 17, 2022

How are the values for HS.Registry.Config set?

In the class HS.Util.Zip.Operations the OnInit method does a get value for **ZipUtility**, **ZipCommand**, **UnZipCommand**.

What is not clear in documentation is where/how those key/values are set. Any assistance would be appreciated.

Method OnInit() As %Status

```
{
    If (##class(HS.Registry.Config).GetKeyValue("/ZipUtility/ZipCommand")="" ||
    (##class(HS.Registry.Config).GetKeyValue("/ZipUtility/UnZipCommand")="" ) {
        $$$LOGERROR($$$Text("Zip Configuration is not configured, please define config
registry entries /ZipUtility/ZipCommand and /ZipUtility/UnZipCommand"))
    }
    Quit $$$OK
}
```

[#Key Value #HealthShare](#)

Product version: IRIS 2020.1

\$ZV: IRIS for Windows (x86-64) 2020.1 (Build 215U) Mon Mar 30 2020 20:14:33 EDT [HealthConnect:2.1.0]

Source URL: <https://community.intersystems.com/post/how-are-values-hsregistryconfig-set>