

Question

[Joseph Tsang](#) · Jun 27, 2022

Populate object by %LoadFromMemory()

Has anyone here been successfully populating a %Persistent object from calling %LoadFromMemory() with a temp. global?

I have tried using Sample.Person in the SAMPLES namespace. It worked EXCEPT that I have no way to set the id into that Sample.Person objec.

```
SAMPLES>Merge array(1) = ^Sample.PersonD(1)
```

```
SAMPLES>Set obj = ##class(Sample.Person).%LoadFromMemory(1,.array)
```

```
SAMPLES>zw obj
```

```
obj=1@Sample.Person ; <OREF>
```

```
+----- general information -----
```

```
|      oref value: 1
```

```
|      class name: Sample.Person
```

```
|      reference count: 2
```

```
+----- attribute values -----
```

```
|      %Concurrency = 1 <Set>
```

```
|      DOB = 32008
```

```
|      Name = "Jaynes,Will W."
```

```
|      SSN = "996-40-3910"
```

```
+----- swizzled references -----
```

```
|      i%FavoriteColors = "" <Set>
```

```
|      r%FavoriteColors = "" <Set>
```

```
|      i%Home = $lb("4693 Maple Drive","Albany","SD",10333) <Set>
```

```
|      r%Home = "" <Set>
```

```
|      i%Office = $lb("3989 First Drive","Miami","KS",76548) <Set>
```

```
|      r%Office = "" <Set>
```

```
|      i%Spouse = ""
```

```
|      r%Spouse = ""
```

```
+----- calculated references -----
```

```
|      Age <Get>
```

```
+-----
```

```
SAMPLES>Write obj.%Id()
```

```
SAMPLES>
```

```
SAMPLES>Do obj.%IdSet(1)
```

```
DO obj.%IdSet(1)
```

```
^
```

```
<PRIVATE METHOD>
```

```
SAMPLES>Write obj.%Id()
```

```
SAMPLES>
```

Populate object by %LoadFromMemory()

Published on InterSystems Developer Community (<https://community.intersystems.com>)

[#Object Data Model](#) [#ObjectScript](#) [#InterSystems IRIS](#)

Product version: IRIS 2020.1

Source URL: <https://community.intersystems.com/post/populate-object-loadfrommemory>