
Question

[Evgeny Shvarov](#) · Oct 6, 2022

What is the Analog of \$DATA for Embedded Python?

Hi folks!

I'm working with a global via Embedded Python as a class method. I init the handler for a global via:

```
gl=iris.gref( "^Global" )
```

What is the way to check the value at index ^Global("x","y")? E.g.:

In ObjectScript I'd do the following:

```
set result=$data(^Global("x","y"))
```

How do I do the same in Embedded Python?

I checked the [documentation](#), but haven't found an answer.

[#Embedded Python](#) [#Globals](#) [#ObjectScript](#) [#InterSystems IRIS](#)

Product version: IRIS 2022.1

Source URL: <https://community.intersystems.com/post/what-analog-data-embedded-python>