
Question

[George Hodder](#) · Oct 14, 2022

Is there a way to flush the routine buffer for a cached routine?

I am having an issue where a routine that is changed is at times still running the old version of code. The routine is compiled after the change.

I am wondering if it is possible that the routine code is still in the routine buffer cache. Is there a way to force a flush or reload of the buffer for a given routine?

[#InterSystems IRIS](#)

Product version: IRIS 2020.1

\$ZV: IRIS for UNIX (Red Hat Enterprise Linux for x86-64) 2020.1 (Build 217120000U) Wed Jun 10 2020

21:40:24 EDT

Source URL: <https://community.intersystems.com/post/there-way-flush-routine-buffer-cached-routine>