

---

Question

[Phillip Wu](#) · Dec 20, 2022

## TrakCare iprint

Hi,

I'm trying to configure iprint for trakcare.

My trakcare runs on linux IRIS(SUSE).

iprint is from MicroFocus [previously Novell]

iprint is a print facility that many printers support where there is a virtual queue for each AD user.

You can log onto any iprint printer using your AD credentials and you can access your virtual print queue. You can then physical print from that virtual queue.

Has anyone configured this for their trakcare?

In order to do this I need to:

- map each trakcare user to their equivalent AD user
- get an iprint driver

However there is no iprint driver for my Linux IRIS.

Has anyone done something similar?

[#TrakCare](#)

Product version: IRIS 2021.2

---

Source URL: <https://community.intersystems.com/post/trakcare-iprint>