

IrisStat options

a	All
b	bits
c	counters
d	dump processes
e	error log
f	global module flags
g	^GLOSTAT info
h	irisstat usage
j	journal
k	prefetch daemons
l	LRU global buffers
m	GFILETAB
n	network
o	clear irisstat
p	processes
q	hibernation semaphores
s	irisstat exe directory
t	run irisstat in loop
u	locks
v	Check versions
w	write daemon queues
B	Blocks in GBFSPECQ
C	inter-job comms
D	sample block collisions
E	cluster status
G	global buffers BDB
H	global buffers SFN/BLK
I	incremental backups
L	license
M	mailbox log
N	ECP
R	routine buffers
S	hang info
T	in-memory tables
V	process memory variables
W	perform thaw
X	device translation table

<https://docs.intersystems.com/iris20191/csp/docbook/DocBook.UI.Page.cls?>

KEY=GCM_irisstat#GCM_irisstat_running_options